

WEBER IRIS KNIFE STRAIGHT, 13CM

Ref. IC60500



26/08/2024 23:02:25 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

WEBER IRIS KNIFE STRAIGHT, 13CM



Ref. IC60500

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040012766	HS Code	90189090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	0.6
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	12.9
Packed unit / Weight	0.05	Bare Unit / Weight (kg)	0.05
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

26/08/2024 23:02:25 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr