

SPO2 SENSOR ADULT CLIP FOR COMEN NC3/NC5/C30

Ref. DBCADUL01



CE
0537

DISPOSITIF MEDICAL
CLASSE
IIb
MEDICAL DEVICE

GARANTIE
1
1 YEAR
GUARANTEE

25/08/2024 21:06:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SPO2 SENSOR ADULT CLIP FOR COMEN NC3/NC5/C30

Ref. DBCADUL01

Designed for accurate and fast measurement of oxygen saturation

Description

Technical specifications: % oxygen: 90-100% accuracy +/-
-2% 70-90%

accuracy +/-3%
<70% not specified

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3661040219417	HS Code	90181990
Warranty (months)	12	Country	Uruguay
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	3
Packed unit / Length (cm)	13	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 25/08/2024 21:06:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr