

X-RAY VIEWER 3 PANELS, 3X45W

Réf. OTA1003



X-RAY VIEWER 3 PANELS, 3X45W



Réf. OTA1003

These negatoscopes feature an X-ray mounting system for easy installation and removal, without affecting the X-ray. They can be placed on a table or mounted on a wall.

- The screen structure is made of 10 cm thick, high-strength aluminum alloy.
- The corners are protected by an ABS-coated steel plate.
- Methacrylic polymer screen
- Negotioscope is in high-strength aluminum alloy
- Interior and exterior coated with antistatic spray powder
- Voltage / Frequency: 220V / 50-60Hz
- Luminance: 2,000 cd/m²
- Average illumination > 6300 lux
- Lighting life: 50,000 h
- Color temperature: 10,000 k
- External dimensions: 116 x 56 x 10 cm
- Screen dimensions: 108 x 43.5 cm

FEATURES			
Medical device	No		
Marquage	CE	Barcode	3661040195872
Guarantee (month):	12	Customs Code	90105010
Sterile	No		
INFORMATIONS		Main unit	Packaged unit
Width (cm) :		Width packaging (cm):	68,5
Length (cm) :		Length packaging (cm) :	122,0
Height (cm) :		Height packaging (cm) :	15,2
Weight (kg) :		Weight packaging (kg) :	16,830
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



11/07/2023 16:53 - This article is available on our website : www.gsh-med.fr with the following reference OTA1003

All characteristics are given for information only and may be changed