

# EXTRA-FLAT X-RAY VIEWER - 1 PANEL

Réf. OTA1101



# EXTRA-FLAT X-RAY VIEWER - 1 PANEL



Réf. OTA1101

These negatoscopes feature an X-ray mounting system for easy installation and removal, with no consequences for the X-ray.

- The screen structure is made of 2.4 cm thick, high-strength aluminum alloy.
- The corners are protected by an ABS-coated steel plate
- Methacrylic polymer screen
- Negotoscope is in high-strength aluminum alloy
- Interior and exterior coated with antistatic spray powder
- Voltage / Frequency: 220V / 30-100Hz
- Average brightness > 1500 cd/m<sup>2</sup> ; >5,000 lux.
- Lighting life: 100,000 h
- Color temperature: 6500 k
- Screen dimensions: 36 x 43.5 cm
- External dimensions: 47 x 54.5 x 2.4 cm

FEATURES			
Medical device	No		
Marquage	CE	Barcode	3661040195889
Guarantee (month):	12	Customs Code	90105010
Sterile	No		
INFORMATIONS		Main unit	Packaged unit
Width (cm) :		Width packaging (cm):	57,0
Length (cm) :		Length packaging (cm) :	67,0
Height (cm) :		Height packaging (cm) :	10,0
Weight (kg) :		Weight packaging (kg) :	6,250
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



11/07/2023 17:03 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference OTA1101

All characteristics are given for information only and may be changed