

TOCO PROBE / WATER PROOF

Ref. BCTOCO15

COMEN



CE
0482



 30/10/2024 22:18:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

TOCO PROBE / WATER PROOF

COMEN

Ref. BCTOCO15

TOCO probe of the cardiotocograph to monitor the contractions of the mother-to-be

Description

IP67 waterproof

Features

CE Marking	CE0482	Medical device	Class IIa
Bar code	3661040219264	HS Code	9018193010
Warranty (months)	6	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/10/2024 22:18:14 - For information purposes only, can be subject to change