

## HARTMANN EAR POLYPUS FORCEPS 8 CM

Ref. IP14508



25/08/2024 20:15:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

# HARTMANN EAR POLYPUS FORCEPS 8 CM



Ref. IP14508

Facilitates the grasping and extraction of foreign bodies in the nose and ear

## Description

Size : 12 x 3 mm

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

100% IN-STORE reference: this reference of service instrumentation is guaranteed to be available immediately.

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040022857	HS Code	90189090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	23	Bare Unit / Width (cm)	8.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.1
Packed unit / Length (cm)	9	Bare Unit / Length (cm)	8
Packed unit / Weight	0.03	Bare Unit / Weight (kg)	0.028
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

25/08/2024 20:15:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)