

LERICHE FORCEPS STRAIGHT NO TEETH 15 CM

Ref. IP17215



25/08/2024 20:17:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

LERICHE FORCEPS STRAIGHT NO TEETH 15 CM



Ref. IP17215

Allows you to perform certain acts in surgery (holding suture threads, compresses, etc.)

Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

100% IN-STORE reference: this reference of service instrumentation is guaranteed to be available immediately.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040023144	HS Code	90189090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	23	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	9	Bare Unit / Length (cm)	15
Packed unit / Weight	0.05	Bare Unit / Weight (kg)	0.035
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

25/08/2024 20:17:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr