

CONTAINER FLOVAC 2L ASPIRET, ASKIR 20, 30, C30BR 31841



Ref. BCASPI103



CE
0123

DISPOSITIF MEDICAL
CLASSE
I_m
MEDICAL DEVICE

GARANTIE
2
2 YEARS
GUARANTEE

13/12/2024 22:12:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

CONTAINER FLOVAC 2L ASPIRET, ASKIR 20, 30, C30BR 31841



Ref. BCASPI103

Description

- FLOVAC® rigid polycarbonate jar / autoclavable container.
- Intended to hold FLOVAC disposable pouches- capacity: 2000ml FLOVAC®

® systems therefore replace the use of autoclavable jars, which reduces cleaning and sterilization costs.

Features

CE Marking	CE0123	Medical device	Class Im (measuring function)
Bar code	3661040161365	HS Code	90189084
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.65
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

13/12/2024 22:12:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr