

REFLEX VERNON HAMMER CHILD RING



Ref. ID06611



 26/08/2024 23:18:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

REFLEX VERNON HAMMER CHILD RING



Ref. ID06611

Description

Spare black ring for reflex hammer from Vernon Child

Total diameter of the ring: 3.73 cm

Rubber width: 0.95 cm

Features

CE Marking	NONE	Medical device	No
Bar code	3661040013596	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	6.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	24.5
Packed unit / Weight	0.01	Bare Unit / Weight (kg)	0.01

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

26/08/2024 23:18:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr