

REFLEX BABINSKY HAMMER ADULT RING



Ref. ID06612



26/08/2024 23:19:40 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

REFLEX BABINSKY HAMMER ADULT
RING

Ref. ID06612

Description

Black Ring

Total diameter of the ring: 4.45 cm

Rubber width: 1.05 cm

Features

CE Marking	NONE	Medical device	No
Bar code	3661040013602	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	4.4
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	12
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	4.4
Packed unit / Weight	0.016	Bare Unit / Weight (kg)	0.014

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

26/08/2024 23:19:40 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr