

DOYEN MOUTH-GAG CHILD, 12 CM

Ref. ID03600



 10/01/2025 00:41:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

DOYEN MOUTH-GAG CHILD, 12 CM



Ref. ID03600

Description

Doyen children's mouth opener in stainless steel with a total length of 12 cm (Hooks: 45 mm - Ends: 21 mm).

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040013442	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	6.3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	3
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	12
Packed unit / Weight	0.081	Bare Unit / Weight (kg)	0.081
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

10/01/2025 00:41:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr