

TUNING FORK WITH CURSOR N°256 VD, 16 CM

Ref. ID26256



10/01/2025 03:06:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

TUNING FORK WITH CURSOR N°256 VD, 16 CM



Ref. ID26256

Makes it possible to carry out a first diagnosis of deafness and to determine its type

Description

Holtex instruments are all made of premium stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to ensure good durability and corrosion resistance.

The steels come into contact with human tissues, so they must comply with international standards and regulations.


Holtex guarantees you medical devices that comply with Directive 93/42/EC and ISO 13485 (2016).


- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" leaflet

Features

| | | | |
|------------------------------------|---------------|-------------------------|-----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3661040013923 | HS Code | 9018.9090 |
| Warranty (months) | 24 | Country | Pakistan |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Bare Unit / Width (cm) | 6 |
| Packed unit / Height (cm) | nc | Bare Unit / Height (cm) | 1.9 |
| Packed unit / Length (cm) | nc | Bare Unit / Length (cm) | 16 |
| Packed unit / Weight | 0.063 | Bare Unit / Weight (kg) | 0.063 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

 10/01/2025 03:06:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr