

SINGLE-ENDED AMALGAM CARRIER DIA. 2.5 MM

Ref. IDPASI25



27/08/2024 00:29:29 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

SINGLE-ENDED AMALGAM CARRIER DIA. 2.5 MM



Ref. IDPASI25

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040015811	HS Code	9018.9090
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	2	Bare Unit / Width (cm)	2
Packed unit / Height (cm)	1	Bare Unit / Height (cm)	1
Packed unit / Length (cm)	15	Bare Unit / Length (cm)	15.5
Packed unit / Weight	0.051	Bare Unit / Weight (kg)	0.051
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 27/08/2024 00:29:29 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr