

BAG FOR OXYGEN CYLINDER FROM 3 TO 5L, RED EB02.016

 **ELITE BAGS**

Ref. METUBE5L_R



 25/08/2024 20:49:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

BAG FOR OXYGEN CYLINDER FROM 3 TO 5L, RED EB02.016



Ref. METUBE5L_R

Description

- Cylindrical bag for oxygen therapy
- Sold empty
- For oxygen bottle of 3 to 5 liters
- Measures: 66 x 35 x 22 Ø cm; weight 1.11 Kg
- Capacity: 30 L
- Material: 1000D Polyester. Washable.

Features

CE Marking	NONE	Medical device	No
Bar code	3661040202938	HS Code	42029270
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	35	Bare Unit / Width (cm)	22
Packed unit / Height (cm)	22	Bare Unit / Height (cm)	35
Packed unit / Length (cm)	66	Bare Unit / Length (cm)	66
Packed unit / Weight	1.11	Bare Unit / Weight (kg)	1.11
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 25/08/2024 20:49:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr