

PROBE COVER PC7200 RADIANT, BOX OF 200 PCS

Ref. CRAPC7200



30/08/2024 16:00:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

PROBE COVER PC7200 RADIANT, BOX OF 200 PCS

Ref. CRAPC7200

Description

Thermometer accessory, for use only with an ear thermometer.

Disposable probe cover, do not use it more than once to avoid cross-contamination.

Compatible with the 4VS Spengler Vital Signs Monitor and ear thermometers such as: Omron Gentle Temp 520 and 521, EDAN M3A...

Can be used with or without dispenser (ref. CRAPC7200 sold without distributor, packaging in French).

Features

CE Marking	CE	Medical device	Class I
Bar code	4710810451546	HS Code	9025900090
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	8	Packed unit / Height (cm)	3
Packed unit / Length (cm)	17	Packed unit / Weight	0.1
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	50
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 16:00:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr