

**POWER SUPPLY (CH,EURO,UK) 12V,
500MA, INPUT 100V-240V (FROM 2017)**

KERN[®]

Ref. EKTRAN01

 30/08/2024 15:39:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

POWER SUPPLY (CH,EURO,UK) 12V, 500MA, INPUT 100V-240V (FROM 2017)



Ref. EKTRAN01

Description

External AC adapter, 100-240 V, 50/60 Hz, Standard EU, CH, UK
Input voltage: 100 V - 240 V, 50/60 Hz

Features

CE Marking	NONE	Medical device	No
Bar code	3661041009925	HS Code	85044090
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	8
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	6
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	10
Packed unit / Weight	nc		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 15:39:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr