

REDON FORCEPS/DISSECTOR 90° 14 CM BABY

Ref. IP49015



29/08/2024 15:37:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

REDON FORCEPS/DISSECTOR 90° 14 CM BABY



Ref. IP49015

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040150918	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	7	Packed unit / Height (cm)	1
Packed unit / Length (cm)	15	Packed unit / Weight	0.038
Bare Unit / Weight (kg)	0.038		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 15:37:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr