

TRANSPORT BAG FOR 4VS

Ref. MDSACT01



30/08/2024 15:54:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

TRANSPORT BAG FOR 4VS



Ref. MDSACT01

Description

- Monitor carrying bag with double zippered side pockets with mesh pocket for storing accessories and a back pocket.
- Transparent crystal face shutter for good visibility of the screen and increases its protection during its transport and use in intervention.
- The monitor is usable without taking it out of its cover, presence of pvc fabric with holes on the top of the bag to allow monitor ventilation.
- All connectors are accessible through an opening in the inner pockets on each side of the device and at the back of the bag.
- Protection of the device guaranteed with a reinforced bottom in black microfiber fabric and with 4 feet of bottom.
- Transport of the bag with shoulder strap.
- Material: Cordura (strong and breathable synthetic fabric)

Features

CE Marking	NONE	Medical device	No
Bar code	3661040218656	HS Code	42029270
Warranty (months)	12	Country	Spain
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	27
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	15
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	42
Packed unit / Weight	nc	Bare Unit / Weight (kg)	0.5
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Bleaching	Bleaching prohibited	Washing	Wash at 60°C
Cleaning	Dry cleaning	Packed unit / Quantity	1

30/08/2024 15:54:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr