


SPO2 SENSOR PEDIATRIC SILICON FOR MASTER PALM 2



Ref. 222255



CE
0537



30/08/2024 16:47:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

SPO2 SENSOR PEDIATRIC SILICON FOR MASTER PALM 2



Ref. 222255

Description

Pediatric silicon spo2 sensor for oximeter MASTER PALM 2

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446030254	HS Code	90181990
Warranty (months)	12	Country	Uruguay
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.084
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	1	Max temperature	50
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 16:47:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr