

## S/S BOWL WITHOUT SPOUT, FLAT BOTTOM, DIAM 100 MM, 300 ML, WITH COVER

Ref. AI99101



 **Holtex**


CE



 29/08/2024 16:15:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

# S/S BOWL WITHOUT SPOUT, FLAT BOTTOM, DIAM 100 MM, 300 ML, WITH COVER



Ref. AI99101

Allows surgical instruments to be stored and transported

## Description

Holtex instruments are all made of premium stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to ensure good durability and corrosion resistance.

The steels come into contact with human tissues, so they must comply with international standards and regulations.

Holtex guarantees you medical devices that comply with Directive 93/42/EC and ISO 13485 (2016).

- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" leaflet

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040160573	HS Code	90189090
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	21	Bare Unit / Width (cm)	21.5
Packed unit / Height (cm)	4	Bare Unit / Height (cm)	4.6
Packed unit / Length (cm)	38	Bare Unit / Length (cm)	38
Packed unit / Weight	0.1	Bare Unit / Weight (kg)	0.096
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 16:15:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)