

SUCTION UNIT CA-MI NEW HOSPIVAC 350, AP FULL 5, 2X5 LITERS AUTOCLAV



Ref. BCASPI507



13/12/2024 22:15:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SUCTION UNIT CA-MI NEW HOSPIVAC 350, AP FULL 5, 2X5 LITERS AUTOCLAV



Ref. BCASPI507

Designed to drain the patient's body fluids, gases, tissues and bones, during or after a procedure. The state-of-the-art oil-free pump system guarantees excellent suction capabilities and maximum vacuum pressure in seconds, without any need for maintenance.

Description

-with foot control and automatic suction inverter

Designed to drain the patient's body fluids, gases, tissues and bones, during or after a procedure. The state-of-the-art oil-free pump system guarantees excellent suction capabilities and maximum vacuum pressure in seconds, without any need for maintenance.

- The precise dashboard and the full range of accessories that come with it make it the ideal device for surgical aspiration and liposuction.
- Accessories included:- Autoclavable

jars with overflow valve

-8x14 Autoclavable Silicone Tubes + Conical connector o 10-11-12mm

-Antibacterial and hydrophobic filter

-Safety container 220 ml

- Technical characteristics

- Engine: Oil-free piston pump requiring no maintenance

-

Power supply: 230V-50Hz

-

Maximum vacuum (adjustable): -0.90 Bar -90 kPa -675 mmHg

- Max free air flow: 60 l/min

- Noise level: 52 dBA

- Weight: 13 kg

- Dimensions: 46 x 85 x 42 cm

 13/12/2024 22:15:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SUCTION UNIT CA-MI NEW HOSPIVAC 350, AP FULL 5, 2X5 LITERS AUTOCLAV



Ref. BCASPI507

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661041009949	HS Code	90189084
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	13
Bare Unit / Weight (kg)	13		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

13/12/2024 22:15:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr