

JARS POLI 2L FOR ASPIRET, ASKIR 20,30 AND C30BR



Ref. BCASPI81



13/12/2024 21:53:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

JARS POLI 2L FOR ASPIRET, ASKIR 20,30 AND C30BR



Ref. BCASPI81

Features

CE Marking	NONE	Medical device	No
Bar code	3661040160979	HS Code	90330090
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.2
Bare Unit / Weight (kg)	0.2		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	1	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

13/12/2024 21:53:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr