

HOLTEX FURNITURE KIT: COUCH+STOOL+STEP STOOL WITH 2 STEPS, BLACK



Ref. OCABREV2_N



29/08/2024 16:46:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

HOLTEX FURNITURE KIT: COUCH+STOOL+STEP STOOL WITH 2 STEPS, BLACK



Ref. OCABREV2_N

Cabinet COMPLETE BLACK - step 2 steps

Description

Set of furniture consisting of:

* 1 stainless steel sofa (OTA003_N)

- With quick disassembly system and adjustable headrest with high resistance automatic rack
- Head and seat padded with very high density foam (30 kg/m3), whose thickness is 5 cm, covered with fireproof washable plasticized fabric (fire resistance classification: M1)
- Removable roller holder (for rolls 50 cm wide) positionable at the head or feet
- maximum load: 150 kg

* 1 stool (OTA007_N)

- Seat ø: 43cm, ep.: 6.5 cm, padded in very high density foam (30kg/m3) with screw-down steel plate.
- Black plastic base with 5 branches ø: 58 cm, very stable with swivel casters to clip-on
- Height adjustable by simple pressure and gas piston
- Height variation from 53 to 70 cm- Maximum load: 150 kg

* 1 step 2 steps (OS0050)

- Structure in stainless steel tube
- Steps covered with non-slip plastic material
- Maximum load: 150 kg

Stainless steel or Epoxy: what's the difference?

Epoxy: the economic choice. The steel is simply coated with a resin to delay corrosion.

Stainless steel: it is steel modified to become completely STAINLESS steel, it does not rust.


And did you know?

Today, more than 90% of stainless steel is made from recycled steel, it's good for the environment!

 29/08/2024 16:46:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

HOLTEX FURNITURE KIT: COUCH+STOOL+STEP STOOL WITH 2 STEPS, BLACK



Ref. OCABREV2_N

Features

CE Marking	NONE	Medical device	No
Bar code	3661040003986	HS Code	9402900000
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	45.6
Bare Unit / Weight (kg)	45.1		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 16:46:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr