

SET OF SILICONE TUBES Ø 6X10 JAR + CONNECTOR FOR ASPIRET, ASKIR 20/30



Ref. BCASPI113

 29/08/2024 15:52:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SET OF SILICONE TUBES Ø 6X10 JAR + CONNECTOR FOR ASPIRET, ASKIR 20/30



Ref. BCASPI113

Features

CE Marking	NONE	Medical device	No
Bar code	3661040151052	HS Code	90330090
Country	Italy		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.05
Bare Unit / Weight (kg)	0.05		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 15:52:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr