

JARS IN POLICARBONATE (MAX 121°C) 400 ML WITH COVER FOR EMIVAC



Ref. BCASPI77



 13/12/2024 22:07:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr

JARS IN POLICARBONATE (MAX 121°C) 400 ML WITH COVER FOR EMIVAC



Ref. BCASPI77

Features

CE Marking	NONE	Medical device	No
Bar code	3661040151076	HS Code	4202121900
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	18.4
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	9.3
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	27
Packed unit / Weight	0.15	Bare Unit / Weight (kg)	0.1
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

13/12/2024 22:07:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr