

WEIGHING PLATFORM WITH 2 BARRIERS (1 FOLDING BARRIER AND DISPLAY HOLDER AND 1 FIXED) MWA (300KG) CLASS III VERIFICATION

KERN®

Ref. EKMWA03_CIII

CE
0297



 27/01/2025 21:50:15 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

WEIGHING PLATFORM WITH 2 BARRIERS (1 FOLDING BARRIER AND DISPLAY HOLDER AND 1 FIXED) MWA (300KG) CLASS III VERIFICATION



Ref. EKMWA03_CIII

Wheelchair weighing platform with 2 barriers (1 folding) and display door approved

Description

Weighing Platform wheelchair

- LCD display, digit height 25 mm
- Display Size (WxDxH) 200 x 55 x 128 mm
- Overall dimensions (LxWxH) 1180 x 1000 x 1130 mm
- Cable length of display approx. 180 cm
- Optional battery operation, 6 x 1.5V AA not included, operating time approx. 20 pm
- Power adapter (external) included • Net weight approx. 57.6 kg

Features

CE Marking	CE0297	Medical device	Class Im (measuring function)
Bar code	3661040173665	HS Code	84231090
Warranty (months)	36	Country	Germany

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	100
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	113
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	118
Packed unit / Weight	94	Bare Unit / Weight (kg)	80

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

27/01/2025 21:50:15 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr