

## BABY SCALE KERN MBC, 20KG

Ref. EKMBC02

# KERN®



CE  
0297


DISPOSITIF MEDICAL  
CLASSE  
**Im**  
MEDICAL DEVICE

GARANTIE  
2 ANS  
2 YEARS  
GUARANTEE

 27/08/2024 08:52:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## BABY SCALE KERN MBC, 20KG



Ref. EKMB02

### Description

- Appreciated in paediatrics thanks to its large ergonomic and removable weighing surface, ideal for infants
- Special function for neonatal services: difference function for calculating weight gain before and after food intake
- Resistant ABS plastic. Hygienic and easy to clean
- Stable and non-slip thanks to its rubber feet
- Hold function: for restless babies lying down, the average value looks for a stable weighing value and "freezes" it
- LCD display with backlight, digit height 25 mm
- Graduation 0,01 kg
- Battery operated (6 x 1.5 AA not included), AC (power supply supplied) or optional with rechargeable battery pack
- Approval possible in Class III
- Medical approval according to 93/42/EEC

Technical specifications:

- Dimensions weighing area :  
L×P×H 600×280×55 mm
- Total dimensions :  
L×P×H 600×400×110 mm
- Net weight approx. 5.5 kg

For any measuring device requiring medical Class III approval, accessories must be ordered at the same time as the device in order for the approval to be valid.

### Features

|                                    |               |                           |                               |
|------------------------------------|---------------|---------------------------|-------------------------------|
| CE Marking                         | CE0297        | Medical device            | Class Im (measuring function) |
| Bar code                           | 4045761095854 | HS Code                   | 84231090                      |
| Warranty (months)                  | 24            | Country                   | Germany                       |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |                               |
| Packed unit / Width (cm)           | nc            | Bare Unit / Width (cm)    | 28                            |
| Packed unit / Height (cm)          | nc            | Packed unit / Length (cm) | nc                            |
| Bare Unit / Length (cm)            | 60            | Packed unit / Weight      | 6                             |
| Bare Unit / Weight (kg)            | 4.5           |                           |                               |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |                               |
| Minimum temperature                | NA            | Max temperature           | NA                            |
| Minimum humidity                   | NA            | Maximum humidity          | NA                            |
| <b>OTHER FEATURES</b>              |               |                           |                               |
| Packed unit / Quantity             | 1             |                           |                               |

27/08/2024 08:52:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)