

KERN MPC-DISPLAY SCALE-250 KG - 100G-CLASS III MED SCALE

Ref. EKMPC01-CIII

KERN®

CE
0297

DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE


GARANTIE
2
2 YEARS
GUARANTEE



 29/08/2024 11:56:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

KERN MPC-DISPLAY SCALE-250 KG - 100G-CLASS III MED SCALE



Ref. EK MPC01-CIII

Bathroom scales EK MPC01 approved in medical class III

Description

- Hygienic and easy to clean
- Stable and non-slip thanks to its height-adjustable rubber feet
- Very stable handle for transporting the scale: Large surface with non-slip and wear-resistant coating for maximum stability
- BMI function for calculating underweight/normal/overweight
- Hold function: for agitated patients, the average value looks for a stable weighing value and "freezes"
- "Mother and child" function: the mother's weight is taken, then the ladder is reset at the touch of a button. Then mother and baby can be weighed together. The child's weight will be automatically displayed on the screen
- Vibration filtering
- High-resolution playback: the measurement can be increased by one decimal place by simply pressing a key
- Large LCD screen, with backlight, height of digits 25 mm
- Graduation 0.1 kg
- Battery operation (6 x 1.5V AA included) or mains (power supply supplied) or optional with rechargeable battery pack
- Approved in class III
- Medical approval according to 93/42/EEC

Technical specifications:

- Display dimensions :L×
P×H 200×130×60 mm
- Net weight approx. 10 kg

For any measuring device requiring medical Class III approval, accessories must be ordered at the same time as the device in order for the approval to be valid.

29/08/2024 11:56:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

KERN MPC-DISPLAY SCALE-250 KG - 100G-CLASS III MED SCALE






Ref. **EKMPC01-CIII**

Features

CE Marking	CE0297	Medical device	Class Im (measuring function)
Bar code	3661040063270	HS Code	84231090
Warranty (months)	24	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	36
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	8
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	36.5
Packed unit / Weight	10.02	Bare Unit / Weight (kg)	8.4
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		


29/08/2024 11:56:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr