

STAGNARA GOUGE (NON SLIPPING HANDLE) 30 CM X 10 MM

Ref. IG09410



 27/08/2024 01:46:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

STAGNARA GOUGE (NON SLIPPING HANDLE) 30 CM X 10 MM



Ref. IG09410

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040018997	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	8	Bare Unit / Width (cm)	8.5
Packed unit / Height (cm)	1	Bare Unit / Height (cm)	1.5
Packed unit / Length (cm)	35	Bare Unit / Length (cm)	35.7
Packed unit / Weight	0.113	Bare Unit / Weight (kg)	0.113
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 01:46:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr