

KILLIAN NASAL SPECULUM 30 MM WITHOUT SCREW

Ref. IS23130



27/08/2024 06:56:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

KILLIAN NASAL SPECULUM 30 MM WITHOUT SCREW



Ref. IS23130

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040031187	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	8.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.7
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	23.5
Packed unit / Weight	0.102	Bare Unit / Weight (kg)	0.102
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 27/08/2024 06:56:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr