

BLUE TOURNIQUET CLIP

Ref. DT05641



27/08/2024 08:51:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

BLUE TOURNIQUET CLIP



Ref. DT05641

Easy to install and quick to open

Description

- Compression, release and opening of the buckle with one hand,
- Progressive release
- Polyamide band, elastic for better comfort.
- ABS closure,
- Its composition and design ensure a long service life even under heavy use.
- Length 45 cm (stretched up to 73 cm)

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040054780	HS Code	63079098
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	5	Bare Unit / Width (cm)	2.5
Packed unit / Height (cm)	3	Packed unit / Length (cm)	10
Bare Unit / Length (cm)	45	Packed unit / Weight	0.03
Bare Unit / Weight (kg)	0.03		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 08:51:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr