

**POWER SUPPLY (CH,EURO,UK) 12V,
500MA, INPUT 100V-240V (BEFORE
2017)**

KERN®

Ref. EKACCE11



29/08/2024 13:38:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

POWER SUPPLY (CH,EURO,UK) 12V, 500MA, INPUT 100V-240V (BEFORE 2017)



Ref. EKACCE11

Description

Input voltage: 100 V - 240 V, 50 / 60 Hz

Power supply/adapter included: EURO | AUS | CH | UK | US

Features

CE Marking	NONE	Medical device	No
Bar code	3661040126807	HS Code	85044095
Warranty (months)	36	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.4
Bare Unit / Weight (kg)	0.4		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

29/08/2024 13:38:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr