

RECHARGEABLE BATTERY PACK, INTERNAL FOR KERN FOB(215X215 MM)&KERN MCB

KERN®

Ref. EKFOB01



 27/08/2024 09:14:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

RECHARGEABLE BATTERY PACK, INTERNAL FOR KERN FOB(215X215 MM)&KERN MCB



Ref. EKFOB01

Description

- Three-color LED battery load display
- Suitable for KERN MPS scales, KERN MWS wheelchair weighing platform, KERN MTS scales with retaining bar, KERN MCB scale chair, KERN MBB baby scale and KERN FOB table scale

Features

CE Marking	NONE	Medical device	No
Bar code	4045761067134	HS Code	84231090
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	8	Packed unit / Height (cm)	3
Packed unit / Length (cm)	11	Packed unit / Weight	0.2
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 27/08/2024 09:14:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
 13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr