

## INTERFACE CABLE FOR PC FOR MODELS KERN ABT, CGB

**KERN**<sup>®</sup>


Ref. EKMP510



 27/08/2024 08:39:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## INTERFACE CABLE FOR PC FOR MODELS KERN ABT, CGB



Ref. EKMP510

770-926 RS-232 interface cable for connecting an external device

### Features

CE Marking	NONE	Medical device	No
Bar code	4045761014824	HS Code	84231090
Country	Germany		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	2	Packed unit / Height (cm)	3
Packed unit / Length (cm)	6	Packed unit / Weight	80
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

27/08/2024 08:39:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)