

OCHSNER TROCART 1=2 FRENCH GAUGE (4 MM)

Ref. IT11612



27/08/2024 08:29:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

OCHSNER TROCART 1=2 FRENCH GAUGE (4 MM)



Ref. IT11612

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040045580	HS Code	9018908500
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	6	Packed unit / Height (cm)	2
Packed unit / Length (cm)	20	Packed unit / Weight	0.05
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 08:29:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr