

WIRE TIGHTENER FOR KIRSCHNER BOW

Ref. IT21100



27/08/2024 07:11:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

WIRE TIGHTENER FOR KIRSCHNER BOW



Ref. IT21100

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040031811	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	9.7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	2.1
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	18.1
Packed unit / Weight	0.326	Bare Unit / Weight (kg)	0.326
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 07:11:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr