

LERICHE RETRACTOR - 28 CM - WIDTH 45 MM

Ref. IV05525



27/08/2024 07:14:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

LERICHE RETRACTOR - 28 CM - WIDTH 45 MM



Ref. IV05525

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040088150	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	21	Bare Unit / Width (cm)	21.6
Packed unit / Height (cm)	4	Bare Unit / Height (cm)	4.6
Packed unit / Length (cm)	28	Bare Unit / Length (cm)	28.1
Packed unit / Weight	0.84	Bare Unit / Weight (kg)	0.401
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 07:14:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr