

**STAND TO ELEVATE THE DISPLAY
DEVICE HEIGHT 45CM; FOR EOB
(315X305MM)**

KERN[®]

Ref. EKEOB06

 29/08/2024 11:45:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

STAND TO ELEVATE THE DISPLAY DEVICE HEIGHT 45CM; FOR EOB (315X305MM)



Ref. EKEOB06

Description

Height 45cm

Features

CE Marking	NONE	Medical device	No
Bar code	4045761019607	HS Code	84231090
Country	Germany		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		


29/08/2024 11:45:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr