

ADSON RETRACTOR 4X5 PRONGS BLUNT 330MM

Ref. IE11235



27/08/2024 08:14:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

ADSON RETRACTOR 4X5 PRONGS BLUNT 330MM



Ref. IE11235

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040044170	HS Code	9018908500
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	10.8
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	3.6
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	32.2
Packed unit / Weight	0.184		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 08:14:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr