


ALUMINIUM BOX 21X10X3 CM, GREY

Ref. AA21103



26/08/2024 17:56:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

ALUMINIUM BOX 21X10X3 CM, GREY



Ref. AA21103

Convenient for storing surgical instruments

Description

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040000015	HS Code	90184990
Warranty (months)	12	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	10	Bare Unit / Width (cm)	10.5
Packed unit / Height (cm)	3	Bare Unit / Height (cm)	3.7
Packed unit / Length (cm)	22	Bare Unit / Length (cm)	21
Packed unit / Weight	0.189	Bare Unit / Weight (kg)	0.189
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

26/08/2024 17:56:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr