


ALUMINIUM BOX 21X10X5 CM, GREY

Ref. AA21105



26/08/2024 17:56:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

ALUMINIUM BOX 21X10X5 CM, GREY



Ref. AA21105

Convenient for storing surgical instruments

Description

- Hot sterilization or autoclave (134°-18 minutes)
- Do not wash in an ultrasonic cleaner or in a washing machine
- Warning: aluminum is incompatible with some decontamination products

Format 21 x 10 x 5 cm

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040000039	HS Code	90184990
Warranty (months)	12	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	10	Bare Unit / Width (cm)	10.3
Packed unit / Height (cm)	5	Bare Unit / Height (cm)	5.7
Packed unit / Length (cm)	22	Bare Unit / Length (cm)	21
Packed unit / Weight	0.211	Bare Unit / Weight (kg)	0.211

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

26/08/2024 17:56:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr