


ALUMINIUM BOX GREY 18X9X3 CM

Ref. AAT1893



26/08/2024 17:57:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

ALUMINIUM BOX GREY 18X9X3 CM



Ref. AAT1893

Convenient for storing surgical instruments

Description

- Hot sterilization or autoclave (134°-18 minutes)
- Do not wash in an ultrasonic cleaner or in a washing machine
- Warning: aluminum is incompatible with some decontamination products

Format 18 x 9 x 3 cm

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040000268	HS Code	90184990
Warranty (months)	12	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	9	Bare Unit / Width (cm)	9.5
Packed unit / Height (cm)	4	Bare Unit / Height (cm)	3
Packed unit / Length (cm)	18	Bare Unit / Length (cm)	18.4
Packed unit / Weight	0.147	Bare Unit / Weight (kg)	0.147

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

26/08/2024 17:57:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr