

S/S BOX 20 X 10 X 6 CM

Ref. AD20106

 **Holtex**



CE



 26/08/2024 18:03:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

S/S BOX 20 X 10 X 6 CM



Ref. AD20106

Convenient for storing surgical instruments

Description

- Sold with lid to protect instruments
- Decontamination with detergent and disinfectant -
- Sterilizable by autoclave (134°-18 minutes)
- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- Available in different sizes

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040013657	HS Code	90189090
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	11	Bare Unit / Width (cm)	10.7
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	6
Packed unit / Length (cm)	21	Bare Unit / Length (cm)	20.5
Packed unit / Weight	0.451	Bare Unit / Weight (kg)	0.451
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

26/08/2024 18:03:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr