

STERILIZING STAINLESS STEEL BOX 25X12X6 CM

Ref. AD25126



25/11/2024 20:44:42 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

STERILIZING STAINLESS STEEL BOX 25X12X6 CM



Ref. AD25126

Convenient for storing surgical instruments

Description

- Sold with lid to protect instruments
- Decontamination with detergent and disinfectant -
- Sterilizable by autoclave (134°-18 minutes)
- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- Available in different sizes

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040013770	HS Code	90189090
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	14	Bare Unit / Width (cm)	12
Packed unit / Height (cm)	8	Bare Unit / Height (cm)	6
Packed unit / Length (cm)	27	Bare Unit / Length (cm)	25
Packed unit / Weight	0.744	Bare Unit / Weight (kg)	0.726
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

25/11/2024 20:44:42 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr