

COTTON BOX 11X14 CM

Ref. AI01114



26/08/2024 18:07:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

COTTON BOX 11X14 CM



Ref. AI01114

Allows cotton to be stored to protect against dust

Description

- Diameter 110 cm
- Sold with lid for effective and hygienic protection
- Decontamination with soapy water (do not use detergent) -
- Sterilizable in autoclave (134°-18 minutes) or in dry air
- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- 100% IN-STORE reference: this service instrumentation reference is guaranteed to be available immediately

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019208	HS Code	90189090
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	12	Bare Unit / Width (cm)	11
Packed unit / Height (cm)	12	Bare Unit / Height (cm)	14
Packed unit / Length (cm)	17	Bare Unit / Length (cm)	14
Packed unit / Weight	0.362	Bare Unit / Weight (kg)	0.362
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

26/08/2024 18:07:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr