

NON-LUBRICATED DISPOSABLE SHEATHS FOR ELECTRONIC THERMOMETERS (1000)

Ref. CT40252



26/08/2024 18:58:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

NON-LUBRICATED DISPOSABLE SHEATHS FOR ELECTRONIC THERMOMETERS (1000)

Ref. CT40252

Description


- Disposable case, double protection, for hygiene when taking temperature.
- The inner case prevents the thermometer from coming into contact with the mucosa.
- The color part facilitates the opening of the inner case to place the thermometer.
- The thin and resistant film guarantees accurate measurement results.
- The lubricant-free model is used for oral and axilar temperature measurement.
- For single use, the outer case protects nursing staff and worktops from contamination.
- Shelf life of the case: 5 years.

Available in lubricated version.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040069104	HS Code	39269097
Country	Sweden		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	15	Bare Unit / Width (cm)	15.5
Packed unit / Height (cm)	8	Bare Unit / Height (cm)	8
Packed unit / Length (cm)	21	Bare Unit / Length (cm)	21
Packed unit / Weight	0.841	Bare Unit / Weight (kg)	0.841
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 26/08/2024 18:58:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr