

## REDON INTRODUCING NEEDLE STRAIGHT (DIAM. 3.5 MM), 15CM

Ref. IA08410



26/08/2024 19:35:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
 [www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## REDON INTRODUCING NEEDLE STRAIGHT (DIAM. 3.5 MM), 15CM



Ref. IA08410

### Features

|                             |               |                         |           |
|-----------------------------|---------------|-------------------------|-----------|
| CE Marking                  | CE            | Medical device          | Class I   |
| Bar code                    | 3661040006413 | HS Code                 | 9018.9090 |
| Warranty (months)           | 24            | Country                 | Pakistan  |
| DIMENSIONS AND WEIGHT       |               |                         |           |
| Packed unit / Width (cm)    | nc            | Bare Unit / Width (cm)  | 11        |
| Packed unit / Height (cm)   | nc            | Bare Unit / Height (cm) | 0.3       |
| Packed unit / Length (cm)   | nc            | Bare Unit / Length (cm) | 15.8      |
| Packed unit / Weight        | 0.018         | Bare Unit / Weight (kg) | 0.018     |
| SPECIFIC STORAGE CONDITIONS |               |                         |           |
| Minimum temperature         | NA            | Max temperature         | NA        |
| Minimum humidity            | NA            | Maximum humidity        | NA        |
| OTHER FEATURES              |               |                         |           |
| Packed unit / Quantity      | 1             |                         |           |

26/08/2024 19:35:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)