

LONDON COLLEGE TWEEZERS 15 CM

Ref. IP04715



27/08/2024 02:42:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LONDON COLLEGE TWEEZERS 15 CM



Ref. IP04715

London College Tweezers 15 cm

Description

London College Tweezers 15 cm

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040021874	HS Code	90189090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.6
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	15
Packed unit / Weight	0.023	Bare Unit / Weight (kg)	0.023
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 02:42:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr