

VALVE BLUNT 5 PCS/CLOWARD RETRACTOR (16 MM X 40,45,50,55,60 MM)

Ref. IV08701



27/08/2024 08:03:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

VALVE BLUNT 5 PCS/CLOWARD RETRACTOR (16 MM X 40,45,50,55,60 MM)



Ref. IV08701

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040042596	HS Code	9018908500
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	3.3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	11.4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	10.4
Packed unit / Weight	0.11	Bare Unit / Weight (kg)	0.11
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 08:03:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr