

LANGE-HOHMANN BONE LEVER 27 CM - 45 MM X 70 MM

Ref. IL11845



27/08/2024 01:58:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

LANGE-HOHMANN BONE LEVER 27 CM - 45 MM X 70 MM



Ref. IL11845

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040020310	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	4	Bare Unit / Width (cm)	4.4
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.4
Packed unit / Length (cm)	26	Bare Unit / Length (cm)	26.8
Packed unit / Weight	0.174	Bare Unit / Weight (kg)	0.174
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 01:58:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr